

EBP07372 - IMMUNIZING PEPTIDE

SKU: EBP07372



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004622.2; NP_150643.2; NP_059992.3; NP_001018064.1
Immunogen	Peptide with sequence C-RNWKRKNTKSMNFD, from the internal region of the protein sequence according to NP_004622.2; NP_150643.2; NP_059992.3; NP_001018064.1.
Product Comments	The peptide was used in the production of EBP07372 - Goat Anti-Apolipoprotein E receptor 2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-RNWKRKNTKSMNFD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Approx 90kDa band observed in Mouse Brain lysates (calculated MW of 96.2kDa according to NP_001074395.1 (mouse). Recommended concentration: 0.3-1µg/ml.

DOCUMENTS

- [Data Sheet](#)