



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP07372 - IMMUNIZING PEPTIDE

SKU: EBP07372

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 004622.2; NP 150643.2; NP 059992.3; NP 001018064.1 ID

Peptide with sequence C-RNWKRKNTKSMNFD, from the internal region of the protein sequence according to **Immunogen**

NP 004622.2; NP 150643.2; NP 059992.3; NP 001018064.1.

The peptide was used in the production of EB07372 - Goat Anti-Apolipoprotein E receptor 2 Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-RNWKRKNTKSMNFD Sequence

Shipping

Product

Instructions Ambient

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

Western Approx 90kDa band observed in Mouse Brain lysates (calculated MW of 96.2kDa according to NP 001074395.1

Blot (mouse). Recommended concentration: 0.3-1µg/ml.