

GOAT ANTI-ARF4 ANTIBODY

SKU: EB05550



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	ADP-ribosylation factor 4 ARF4
Names	
Accession ID	NP_001651
Immunogen	Peptide with sequence GLTISSLFSRL-C, from the N Terminus of the protein sequence according to NP_001651.
Peptide Sequence	GLTISSLFSRL-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Mouse
Human Gene ID	378
Mouse Gene ID	11843
Rat Gene ID	79120
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 20-22kDa band seen in Mouse Brain lysate. Recommended for use at 1-3µg/ml.