

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ARF4 ANTIBODY

SKU: EB05550

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias ADP-ribosylation factor 4|ARF4

Names

Accession NP_001651

Immunogen Peptide with sequence GLTISSLFSRL-C, from the N Terminus of the protein sequence according to NP_001651.

Peptide

GLTISSLFSRL-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat **Species**

Reactive

Mouse **Species**

Human 378

Gene ID

Mouse 11843

Gene ID Rat Gene ID 79120

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

Western

Western blot: Approx 20-22kDa band seen in Mouse Brain lysate. Recommended for use at 1-3µg/ml. Blot