

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ARHGEF5 ANTIBODY

SKU: EB05331

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

transforming immortalized mammary oncogene|guanine nucleotide regulatory protein TIM|oncogene **Alias**

TIM|TIM1|GEF5|TIM|P60|Rho quanine nucleotide exchange factor (GEF) 5|ARHGEF5 **Names**

Usage

Summary

Accession

NM_005435

Blocking

EBP05331 **Peptide**

Peptide with sequence C-VKTAKLQLVEQQA, from the C Terminus of the protein sequence according to **Immunogen**

NM 005435.

Peptide

C-VKTAKLQLVEQQA Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Human, Mouse **Species**

Reactive

Human **Species**

Human

7984 Gene ID

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Western

Western blot: Approx 60-65kDa band seen in 293 lysate. Recommended for use at 1-3μg/ml. **Blot**