

GOAT ANTI-ARHGEF5 ANTIBODY

SKU: EB05331



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	transforming immortalized mammary oncogene guanine nucleotide regulatory protein TIM oncogene TIM TIM1 GEF5 TIM P60 Rho guanine nucleotide exchange factor (GEF) 5 ARHGEF5
Usage Summary	raju
Accession ID	NM_005435
Blocking Peptide	EBP05331
Immunogen	Peptide with sequence C-VKTAKLQLVEQQA, from the C Terminus of the protein sequence according to NM_005435.
Peptide Sequence	C-VKTAKLQLVEQQA
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	7984
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 60-65kDa band seen in 293 lysate. Recommended for use at 1-3µg/ml.