

GOAT ANTI-ARID3C ANTIBODY

SKU: EB08815



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	ARID3C AT rich interactive domain 3C (BRIGHT-like) AT rich interactive domain 3C (BRIGHT- like)
Accession ID	NP_001017363.1
Blocking Peptide	EBP08815
Immunogen	Peptide with sequence C-REKLAPEEPPEKR, from the internal region of the protein sequence according to NP_001017363.1.
Peptide Sequence	C-REKLAPEEPPEKR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Human Gene ID	138715
Mouse Gene ID	550619
Rat Gene ID	502946
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Preliminary experiments in Human Lymph Node, Spleen, Tonsil and Thymus lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?

Application
Type Pep-ELISA

DOCUMENTS

- [Data Sheet](#)