



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ARL1 ANTIBODY

SKU: EB05542

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias ARFL1|ADP-ribosylation factor-like 1|ARL1

Names

Accession NP_001168 ID

Blocking

EBP05542 **Peptide**

Peptide with sequence C-EAMEWLVETLKSRQ, from the C Terminus of the protein sequence according to

Immunogen NP 001168.

Peptide

C-EAMEWLVETLKSRQ Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Human, Mouse, Rat **Species**

Reactive

Human **Species**

Human

400

Gene ID

Mouse

66941 **Gene ID**

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Western Western blot: Approx 20kDa band seen in Human Pancreas lysate [Predicted MW of 20kDa according to

NP_001168]. Recommended for use at 1-3µg/ml. **Blot**