

GOAT ANTI-ARL1 ANTIBODY

SKU: EB05542



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	ARFL1 ADP-ribosylation factor-like 1 ARL1
Accession ID	NP_001168
Blocking Peptide	EBP05542
Immunogen	Peptide with sequence C-EAMEWLVELKSRQ, from the C Terminus of the protein sequence according to NP_001168.
Peptide Sequence	C-EAMEWLVELKSRQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	400
Mouse Gene ID	66941
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 20kDa band seen in Human Pancreas lysate [Predicted MW of 20kDa according to NP_001168]. Recommended for use at 1-3µg/ml.