

GOAT ANTI-ARL4D / ARL4L ANTIBODY

SKU: EB05547



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	ADP-ribosylation factor-like 4D ARL4D ADP-ribosylation factor-like 6 ARL6 ADP-ribosylation factor 4-like ARF4L
Accession ID	NP_001652.2
Immunogen	Peptide with sequence GNHLTEMAPTASS-C, from the N Terminus of the protein sequence according to NP_001652.2.
Peptide Sequence	GNHLTEMAPTASS-C
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Reactive Species	Human
Human Gene ID	379
Mouse Gene ID	80981
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western Blot: Approx 23-25kDa band observed in Human Lung and Placenta lysates (predicted MW of 23kDa according to NP_001652). Recommended for use at 1-3µg/ml.
Application Type	Pep-ELISA, WB