


GOAT ANTI-AHR (AA749-763) ANTIBODY

SKU: EB11767



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	class E basic helix-loop-helix protein 76 bHLHe76 aryl hydrocarbon receptor aromatic hydrocarbon receptor AH-receptor AHR ah receptor
Accession ID	NP_001612.1
Blocking Peptide	EBP11767
Immunogen	Peptide with sequence C-KHGLNPQSAIITPQT, from the internal region of the protein sequence according to NP_001612.1.
Peptide Sequence	C-KHGLNPQSAIITPQT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	196
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 100kDa band observed in lysates of cell line HeLa (calculated MW of 96.1kDa according to NP_001612.1). Recommended concentration: 0.2-0.5µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

