

GOAT ANTI-ARYLSULFATASE C / STS ANTIBODY

SKU: EB08062



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa

25kDa
20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	STS arylsulfatase C steroid sulfatase (microsomal), isozyme S ARSC ARSC1 ASC ES SSDD estrone sulfatase steryl-sulfatase steryl-sulfate sulfohydrolase
Accession ID	NP_000342.2
Blocking Peptide	EBP08062
Immunogen	Peptide with sequence C-QAGQKIDEPTSN, from the internal region of the protein sequence according to NP_000342.2.
Peptide Sequence	C-QAGQKIDEPTSN
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	412
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 70kDa band observed in Human Placenta lysates (calculated MW of 65.5kDa according to NP_000342.2). Recommended concentration: 0.01-0.03µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

