Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-ASCL1A (ZEBRAFISH) ANTIBODY**

**SKU:** EB07951

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 ua

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / pituitary absent|pia|achaete-scute complex-like

Alias 

**Names** complex-like 1a (Drosophila)|ascl1a

Accession

NP\_571294.1

**Blocking Peptide** 

EBP07951

Peptide with sequence QCSNKSASKQAKRQR, from the internal region of the protein sequence according to NP 571294.1 **Immunogen** 

**Peptide** Sequence

**OCSNKSASKOAKROR** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the

Method

immunizing peptide.

Refrigerated Instructions

**Predicted Species** 

Zebrafish

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

Detection Antibody detection limit dilution 1:8000.

Limit

Western Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product.

Therefore we cannot recommend an optimal concentration and the product is investigative grade.

**Application** 

Type

**Blot** 

Pep-ELISA