

GOAT ANTI-ASCL1A (ZEBRAFISH) ANTIBODY

SKU: EB07951



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	pituitary absent[<i>pia</i>]achaete-scute complex-like 1a zasha zash1a si:dkey-264k15.2 im:7138587 id:ibd5085 asha ash Zash-1a MGC109810 DKEY-264K15.2 achaete-scute complex-like 1a (<i>Drosophila</i>) ascl1a
Accession ID	NP_571294.1
Blocking Peptide	EBP07951
Immunogen	Peptide with sequence QCSNKSASKQAKRQR, from the internal region of the protein sequence according to NP_571294.1
Peptide Sequence	QCSNKSASKQAKRQR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Zebrafish
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

