



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

GOAT ANTI-ASF1A ANTIBODY

SKU: EB08162





SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	anti-silencing function 1A CCG1-interacting factor A ASF1 anti-silencing function 1 homolog A HSPC146 DKFZP547E2110 CGI-98, CIA RP3-329L24.1 ASF1 anti-silencing function 1 homolog A (S. cerevisiae) ASF1A
Accession ID	NP_054753.1
Blocking Peptide	EBP08162
Immunogen	Peptide with sequence C-ENSLNVMLESH, from the C Terminus of the protein sequence according to NP_054753.1.
Peptide Sequence	C-ENSLNVMLESH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	25842
Mouse Gene ID	66403
Rat Gene ID	294408
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 21kDa band observed in lysates of cell line Jurkat (calculated MW of 23.0kDa according to NP_054753.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB



www.everestbiotech.com

Email: customerservice@vectorlabs.com

Telephone: [\(650\) 697-3600](tel:(650)697-3600)

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

