

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ASF1A / HSPC146 ANTIBODY

SKU: EB07788



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

anti-silencing function 1A|CCG1-interacting factor A|ASF1 anti-silencing function 1 homolog Synonym / A|DKFZP547E2110|CIA|CGI-98|RP3-329L24.1|ASF1 anti-silencing function 1 homolog A (S.

cerevisiae)|HSPC146|ASF1A Names

Accession

Alias

ID

NP_054753.1

Blocking

EBP07788 **Peptide**

Peptide with sequence C-EDNTEKLEDAE, from the internal region of the protein sequence according to **Immunogen**

NP 054753.1.

Peptide

C-EDNTEKLEDAE Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog **Species**

Reactive

Human **Species**

Human

25842 Gene ID

Mouse

66403 **Gene ID**

Rat Gene ID 294408

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade ELISA

Detection Antibody detection limit dilution 1:1000.

Limit

Western Approx. 23kDa band observed in lysates of cell line Jurkat (calculated MW of 23.0kDa according to

NP_054753.1). Recommended concentration: 0.01-0.03µg/ml. Blot

Application

Pep-ELISA, WB **Type**

> Goat Anti-ASF1A / HSPC146 Antibody https://everestbiotech.com/products/goat-anti-asfla-hspc146-antibody/





Telephone: (650) 697-3600

GALLERY IMAGES

e erest

250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa