


GOAT ANTI-ASF1A / HSPC146 ANTIBODY

SKU: EB07788



250kDa
150kDa
100kDa
75kDa

50kDa

37kDa

25kDa
20kDa

15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names anti-silencing function 1A|CCG1-interacting factor A|ASF1 anti-silencing function 1 homolog A|DKFZP547E2110|CIA|CGI-98|RP3-329L24.1|ASF1 anti-silencing function 1 homolog A (S. cerevisiae)|HSPC146|ASF1A

Accession ID NP_054753.1

Blocking Peptide EBP07788

Immunogen Peptide with sequence C-EDNTEKLEDAE, from the internal region of the protein sequence according to NP_054753.1.

Peptide Sequence C-EDNTEKLEDAE

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human, Mouse, Rat, Dog

Reactive Species Human

Human Gene ID 25842

Mouse Gene ID 66403

Rat Gene ID 294408

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA Detection Limit Antibody detection limit dilution 1:1000.

Western Blot Approx. 23kDa band observed in lysates of cell line Jurkat (calculated MW of 23.0kDa according to NP_054753.1). Recommended concentration: 0.01-0.03µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

