

GOAT ANTI-ATF1 ANTIBODY

SKU: EB05196



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	cAMP-dependent transcription factor-1 RNA binding protein-activating transcription factor-1 fusion protein TREB36 FUS/ATF-1 EWS-ATF1 activating transcription factor 1 ATF1
Accession ID	NP_005162.1
Blocking Peptide	EBP05196
Immunogen	Peptide with sequence C-ELKTLKDLYSNKS, from the C Terminus of the protein sequence according to NP_005162.1.
Peptide Sequence	C-ELKTLKDLYSNKS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	466
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western blot: Approx 37kDa band seen in A431 cell lysates. Recommended for use at 0.5-2µg/ml.
Application Type	Pep-ELISA, WB

