Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ATF1 ANTIBODY

SKU: EB05196

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

cAMP-dependent transcription factor-1|RNA binding protein-activating transcription factor-1 fusion **Alias** protein|TREB36|FUS/ATF-1|EWS-ATF1|activating transcription factor 1|ATF1

Names

Accession NP_005162.1 ID

Blocking

EBP05196 **Peptide**

Peptide with sequence C-ELKTLKDLYSNKS, from the C Terminus of the protein sequence according to

Immunogen NP 005162.1.

Peptide

C-ELKTLKDLYSNKS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human

Reactive

Human **Species**

Human

466 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:32000.

Limit

Western Western blot: Approx 37kDa band seen in A431 cell lysates. Recommended for use at 0.5-2µg/ml. **Blot**

Application

Pep-ELISA, WB

Type





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600