




Telephone: (650) 697-3600

GOAT ANTI-ATF6 ANTIBODY

SKU: EB05245

e erest



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

ATF6A|activating transcription factor 6|ATF6 **Alias**

Names

Usage

Immunofluorescence: Strong expression of the protein seen in the nuclei and cytoplasm of A549 and in the endoplasmic reticulum and cytoplasm of HeLa cells. Recommended concentration: 10µg/ml.

Flow Cytometry: Flow cytometric analysis of HeLa cells. Recommended concentration: **Summary**

10ug/ml.

Accession

NP_031374.2


Blocking

ID

EBP05245 Peptide











Peptide with sequence C-EATHVVSTIPESLQ, from the C or N Terminus of the protein sequence according to **Immunogen** NP 031374.2.

Peptide C-EATHVVSTIPESLQ **Sequence**

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human

Species Human

22926 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Antibody detection limit dilution 1:8000.

Limit

Application

Pep-ELISA, IF, FC **Type**

GALLERY IMAGES

