

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-ATP11B ANTIBODY

SKU: EB05440

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

DKFZP434N1615|DKFZP434J238|KIAA0956|ATPIR|ATPIF|ATPase, Class VI, type 11B|ATP11B **Alias**

Names

Immunogen Peptide with sequence DRSILTLSTMDSSTC, from the C Terminus of the protein sequence according to .

Peptide

DRSILTLSTMDSSTC Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted Human **Species**

Reactive

Human **Species**

Human 23200

Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite medium.png

Western

Western blot: Approx 85-90kDa band seen in 293 lysate. Recommended for use at 1-3μg/ml. **Blot**