

GOAT ANTI-ATP11B ANTIBODY

SKU: EB05440



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	DKFZP434N1615 DKFZP434J238 KIAA0956 ATPIR ATPIF ATPase, Class VI, type 11B ATP11B
Names	
Immunogen	Peptide with sequence DRSILTLSTMDSSTC, from the C Terminus of the protein sequence according to .
Peptide Sequence	DRSILTLSTMDSSTC
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	23200
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 85-90kDa band seen in 293 lysate. Recommended for use at 1-3µg/ml.