


GOAT ANTI-AVPR1B ANTIBODY

SKU: EB08785



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	vasopressin V1B receptor pituitary vasopressin receptor 3 arginine vasopressin receptor 3 antidiuretic hormone receptor 1B AVPR3 arginine vasopressin receptor 1B AVPR1B
Accession ID	NP_000698.1
Blocking Peptide	EBP08785
Immunogen	Peptide with sequence C-PRDLELADGEGTAE, from the C Terminus of the protein sequence according to NP_000698.1.
Peptide Sequence	C-PRDLELADGEGTAE
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human

Human Gene ID	553
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	In paraffin embedded Human Pituitary shows bright cytoplasm staining of a few active secretory cells in the anterior pituitary. Recommended concentration, 5-10µg/ml.
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx. 170kDa band observed in Human Brain (Substantia Nigra) lysates (calculated MW of 47.0kDa according to NP_000698.1). The observed molecular weight corresponds to findings with EB08786, a product of different design reacting to the Mouse protein. Recommended concentration: 0.05-0.2µg/ml.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

