




Telephone: (650) 697-3600

## **GOAT ANTI-AVPR1B ANTIBODY**

**SKU:** EB08785



## **SPECIFICATIONS**

**Unit Size** 100 μg

Synonym /

Alias Names vasopressin V1B receptor|pituitary vasopressin receptor 3|arginine vasopressin receptor 3|antidiuretic

hormone receptor 1B|AVPR3|arginine vasopressin receptor 1B|AVPR1B

Accession

ID

NP\_000698.1

**Blocking Peptide** 

EBP08785

**Immunogen** 

Peptide with sequence C-PRDLELADGEGTAE, from the C Terminus of the protein sequence according to NP 000698.1.

**Peptide** 

Sequence

C-PRDLELADGEGTAE

Shipping

Instructions

Refrigerated

**Predicted** 

**Species** 

Human

Reactive **Species** 

Human





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

Human Gene ID

553

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**IHC Results** 

In paraffin embedded Human Pituitary shows bright cytoplasm staining of a few active secretory cells in the

anterior pituitary. Recommended concentration, 5-10µg/ml.

**ELISA** 

**Detection** Antibody detection limit dilution 1:32000.

Limit Western

Approx. 170kDa band observed in Human Brain (Substantia Nigra) lysates (calculated MW of 47.0kDa according to NP 000698.1). The observed molecular weight corresponds to findings with EB08786, a product

of different design reacting to the Mouse protein. Recommended concentration: 0.05-0.2µg/ml.


Application

Type

**Blot** 

Pep-ELISA, WB, IHC

## **GALLERY IMAGES**









Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600