

GOAT ANTI-AVPR1B (MOUSE) ANTIBODY

SKU: EB08786



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	VPR3 V3/V1b pituitary vasopressin receptor AVPR3, V3/V1b, VIBR, VPR vasopressin V1B receptor pituitary vasopressin receptor 3 arginine vasopressin receptor 3 antidiuretic hormone receptor 1B AVPR3 arginine
Names	vasopressin receptor 1B AAvpr1b
Accession ID	NP_036054.1
Blocking Peptide	EBP08786
Immunogen	Peptide with sequence C-KPAGSLKDLEQVD, from the internal region (near C Terminus) of the protein sequence according to NP_036054.1.
Peptide Sequence	C-KPAGSLKDLEQVD
Shipping Instructions	Refrigerated
Predicted Species	Mouse, Rat
Reactive Species	Rat
Mouse Gene ID	26361
Rat Gene ID	29462
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx. 170kDa band observed in Rat Brain lysates (calculated MW of 46.5kDa according to NP_036054.1). The observed molecular weight corresponds to findings with EB08785, a product of different design reacting to the C-terminus of the Human protein. Recommended concentration: 0.05-0.2µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

