


GOAT ANTI-BACH1 (AA250-263) ANTIBODY

SKU: EB12081



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	BACH1 BTB and CNC homology 1, basic leucine zipper transcription factor 1 BACH-1 BTB and CNC homolog 1 HA2303 basic region leucine zipper transcriptional regulator BACH1 transcription regulator protein BACH1
Accession ID	NP_001177.1
Blocking Peptide	EBP12081
Immunogen	Peptide with sequence HASVQPNERSENEC, from the internal region of the protein sequence according to NP_001177.1.
Product Comments	Reported variants represent identical protein: NP_996749.1, NP_001177.1
Peptide Sequence	HASVQPNERSENEC
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	571
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 80kDa band observed in Human Spleen and Peripheral Blood Lymphocytes lysates (calculated MW of 82.0kDa according to NP_001177.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

