

GOAT ANTI-BAG3 / BIS/ CAIR1 ANTIBODY

SKU: EB05486



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	BAG-family molecular chaperone regulator-3 BCL2-binding athanogene 3 docking protein CAIR-1 Bcl-2-binding protein DKFZp434E0610 MGC104307 CAIR-1 BAG-3 BIS BCL2-associated athanogene 3 BAG3
Accession ID	NP_004272.2
Blocking Peptide	EBP05486
Immunogen	Peptide with sequence C-SSMTDTPGNPAAP, from the C Terminus of the protein sequence according to NP_004272.2.
Peptide Sequence	C-SSMTDTPGNPAAP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	9531
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Kidney. Recommended concentration: 7µg/ml.
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Approx 80-85kDa band observed in lysates of cell line MCF7 (calculated MW of 61.6kDa according to NP_004272.2). This molecular weight is routinely observed by other sources. Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

