

GOAT ANTI-BAIAP2 ANTIBODY

SKU: EB05176

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

10kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	insulin receptor substrate p53/p58 fas ligand-associated factor 3 brain-specific angiogenesis inhibitor 1-associated protein 2 FLAF3 insulin receptor substrate p53 IRSP53 BAP2 BAI1-associated protein 2 BAIAP2
Accession ID	NP_059345.1
Blocking Peptide	EBP05176
Immunogen	Peptide with sequence C-HGDGSARTLAGR, from the C Terminus of the protein sequence according to NP_059345.1.
Product Comments	This antibody is expected to recognize isoform 2 (NP_059345.1) only.
Peptide Sequence	C-HGDGSARTLAGR
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	10458
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx 60kDa band observed in Human Placenta lysates (calculated MW of 60.9kDa according to NP_059345.1). Recommended concentration: 0.2-0.5µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

