


GOAT ANTI-BAP / SIL1 ANTIBODY

SKU: EB06071



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names SIL1 protein|ULG5|MSS|SIL1 homolog, endoplasmic reticulum chaperone (S. cerevisiae)|BiP-associated protein|endoplasmic reticulum chaperone SIL1, homolog of yeast|BAP|SIL1

Accession ID NP_071909.1; NP_001032722.1

Blocking Peptide EBP06071

Immunogen Peptide with sequence C-ELLGSVNSLLKELR, from the C Terminus of the protein sequence according to NP_071909.1; NP_001032722.1.

Product Comments Variants (NP_071909.1; NP_001032722.1) encode the same isoform.

Peptide Sequence C-ELLGSVNSLLKELR

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 64374

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results Paraffin embedded Human Liver and Testis. Recommended concentration: 5µg/ml.

ELISA Detection Limit Antibody detection limit dilution 1:8000.

Western Blot Approx 55-60 kDa band observed in Human Liver and Kidney extracts (predicted MW of 53kDa according to NP_071909.1). Recommended concentration: 1-2µg/ml.

Application Type Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

