

GOAT ANTI-BAP29 ANTIBODY

SKU: EB05402



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	B-cell receptor-associated protein BAP29 BAP29
Names	
Accession ID	NP_061332
Blocking Peptide	EBP05402
Immunogen	Peptide with sequence CELQDRLERGNKKRL, from the C Terminus of the protein sequence according to NP_061332.
Peptide Sequence	CELQDRLERGNKKRL
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Mouse
Human Gene ID	55973
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 25kDa band seen in Mouse Spleen lysate. Recommended for use at 1-3µg/ml.