

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-BAP29 ANTIBODY

SKU: EB05402

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias B-cell receptor-associated protein BAP29|BAP29

Names

Accession NP_061332 ID

Blocking

EBP05402

Peptide

Peptide with sequence CELQDRLERGNKKRL, from the C Terminus of the protein sequence according to **Immunogen**

NP 061332.

Peptide

CELQDRLERGNKKRL Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Mouse **Species**

Human

55973

Gene ID

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Grade Western

Western blot: Approx 25kDa band seen in Mouse Spleen lysate. Recommended for use at 1-3µg/ml. **Blot**