

GOAT ANTI-BARD1 ANTIBODY

SKU: EB05814



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	BRCA1-associated RING domain gene 1 BRCA1-associated RING domain 1 BRCA1 associated RING domain 1 BARD1
Accession ID	NP_000456
Blocking Peptide	EBP05814
Immunogen	Peptide with sequence CVMSFELLPLDS, from the C Terminus of the protein sequence according to NP_000456.
Product Comments	This antibody is expected to recognise both human isoforms of this protein.
Peptide Sequence	CVMSFELLPLDS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	580
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Jurkat, U937, 293, NCI-H460 and Human Testis lysates (not nuclear preparations) at up to 1µg/ml.