



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-BARD1 ANTIBODY**

**SKU:** EB05814

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

BRCA1-associated RING domain gene 1|BRCA1-associated RING domain 1|BRCA1 associated RING domain Alias 1|BARD1

**Names** 

Accession NP 000456 ID

**Blocking** 

**Peptide** 

EBP05814

Immunogen Peptide with sequence CVMSFELLPLDS, from the C Terminus of the protein sequence according to NP 000456.

**Product** 

This antibody is expected to recognise both human isoforms of this protein. Comments

**Peptide** 

Sequence

**CVMSFELLPLDS** 

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Instructions

Refrigerated

**Predicted Species** 

Human

Human

580

**Gene ID** 

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

**ELISA** 

**Detection** 

Antibody detection limit dilution 1:32000.

Limit

Western Western Blot: No signal obtained yet but low background observed in Jurkat, U937, 293, NCI-H460 and Human

**Blot** Testis lysates (not nuclear preparations) at up to 1µg/ml.