

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-BAT8 / G9A ANTIBODY

SKU: EB05722

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

ankyrin repeat-containing protein|G9A histone methyltransferase|NG36 protein|NG36/G9a|NG36|G9A|HLA-B **Alias**

associated transcript 8|BAT8 **Names**

Accession

NP 006700; NP 079532

Peptide with sequence C-LPELGSLPPVNT, from the C Terminus of the protein sequence according to NP_006700; Immunogen NP_079532.

Product Comments

This antibody is expected to recognise both human isoforms of this protein.

Peptide

C-LPELGSLPPVNT Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

10919 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:2000.

Limit

Western Western Blot: No signal obtained yet but low background observed in Jurkat, HepG2 and U937 lysates at upto

Blot $1\mu g/ml$.