

GOAT ANTI-BAT8 / G9A ANTIBODY

SKU: EB05722



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	ankyrin repeat-containing protein G9A histone methyltransferase NG36 protein NG36/G9a NG36 G9A HLA-B associated transcript 8 BAT8
Accession ID	NP_006700; NP_079532
Immunogen	Peptide with sequence C-LPELGSLPPVNT, from the C Terminus of the protein sequence according to NP_006700; NP_079532.
Product Comments	This antibody is expected to recognise both human isoforms of this protein.
Peptide Sequence	C-LPELGSLPPVNT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	10919
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:2000.
Western Blot	Western Blot: No signal obtained yet but low background observed in Jurkat, HepG2 and U937 lysates at upto 1µg/ml.