

GOAT ANTI-BBC3 ANTIBODY

SKU: EB10389



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	BBC3 PUMA JFY1 Bcl-2 binding component 3 BCL2 binding component 3
Accession ID	NP_001120712.1
Blocking Peptide	EBP10389
Immunogen	Peptide with sequence C-DGGRPQRTVR, from the internal region of the protein sequence according to NP_001120712.1.
Product Comments	This antibody is expected to recognize isoform 1 (NP_001120712.1) only.
Peptide Sequence	C-DGGRPQRTVR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	27113
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 24kDa band observed in Human peripheral blood mononucleocytes lysates (calculated MW of 26.5kDa according to NP_001120712.1). Recommended concentration: 0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

