



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-BBC3 ANTIBODY

SKU: EB10389





Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias BBC3|PUMA|JFY1|Bcl-2 binding component 3|BCL2 binding component 3

Names

Accession

NP_001120712.1

Blocking

EBP10389 **Peptide**

Peptide with sequence C-DGGRPQRTVR, from the internal region of the protein sequence according to **Immunogen**

NP_001120712.1.

Product Comments

This antibody is expected to recognize isoform 1 (NP 001120712.1) only.

Peptide

Sequence

C-DGGRPQRTVR

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human

Species Human

27113 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Western Approx 24kDa band observed in Human peripheral blood mononucleocytes lysates (calculated MW of 26.5kDa

according to NP_001120712.1). Recommended concentration: 0.1-0.3µg/ml. Blot

Application

Pep-ELISA, WB Type

GALLERY IMAGES







250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa