

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-BBS1 (MOUSE) ANTIBODY

SKU: EB10601

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Bardet-Biedl syndrome 1 homolog| Bardet-Biedl syndrome 1| D19Ertd609e| Al451249| Bardet-Biedl syndrome **Alias**

1|Bbs1 **Names**

Accession NP 001028300.1

Blocking

EBP10601 **Peptide**

Peptide with sequence C-NGSIYILRRDSKHPK, from the internal region of the protein sequence according to Immunogen

NP 001028300.1.

Peptide

C-NGSIYILRRDSKHPK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog **Species**

Human

582 Gene ID

Mouse

Gene ID

52028

Rat Gene ID 309156

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

ELISA

Detection Antibody detection limit dilution 1:4000.

Limit

Application Pep-ELISA

Type