

GOAT ANTI-BBS1 (MOUSE) ANTIBODY

SKU: EB10601



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Bardet-Biedl syndrome 1 homolog Bardet-Biedl syndrome 1 D19Ert609e AI451249 Bardet-Biedl syndrome 1 Bbs1
Accession ID	NP_001028300.1
Blocking Peptide	EBP10601
Immunogen	Peptide with sequence C-NGSIYILRRDSKHPK, from the internal region of the protein sequence according to NP_001028300.1.
Peptide Sequence	C-NGSIYILRRDSKHPK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Human Gene ID	582
Mouse Gene ID	52028
Rat Gene ID	309156
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:4000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)