

GOAT ANTI-BBS4 (AA 33-46) ANTIBODY

SKU: EB10604



250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	BBS4 Bardet-Biedl syndrome 4 Bardet-Biedl syndrome 4 protein
Accession ID	NP_149017.2
Blocking Peptide	EBP10604
Immunogen	Peptide with sequence C-EKQNWLIHLHYIRK, from the Internal Region (near N Terminus) of the protein sequence according to NP_149017.2.
Peptide Sequence	C-EKQNWLIHLHYIRK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Dog, Pig
Reactive Species	Human
Human Gene ID	585
Mouse Gene ID	102774
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Approx 65kDa band observed in Human Testis lysates (calculated MW of 58.3kDa according to NP_149017.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

