

GOAT ANTI-BCAP31 / BAP31 ANTIBODY

SKU: EB05397

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	p28 Bap31 OTTHUMP00000025978 BCR-associated protein Bap31 6C6 antigen CDM B-cell receptor-associated protein 31 6C6-AG CDM accessory protein BAP31 DXS1357E BAP31 BCAP31
Accession ID	NP_001132929.1; NP_005736.3
Blocking Peptide	EBP05397
Immunogen	Peptide with sequence C-QAAVDGPMDDKKEE, from the C Terminus of the protein sequence according to NP_001132929.1; NP_005736.3.
Product Comments	This antibody is expected to recognise isoform a (NP_001132929.1) and isoform b (NP_005736.3; NP_001132913.1).
Peptide Sequence	C-QAAVDGPMDDKKEE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Pig, Cow
Reactive Species	Human
Human Gene ID	10134
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Kidney and Brain (Cerebellum and Cortex). Recommended concentration: 3.75µg/ml.
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 30kDa band observed in Human Placenta and Brain lysate (predicted MW of 28kDa according to NP_005736). Recommended for use at 1-3µg/ml.
Application Type	Pep-ELISA, WB, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



