

GOAT ANTI-BCL2L15 ANTIBODY

SKU: EB10383



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	BCL2L15 pro-apoptotic Bcl-2 protein FLJ22588 C1orf178 Bfk BCL2-like 15
Accession ID	NP_001010922.1
Blocking Peptide	EBP10383
Immunogen	Peptide with sequence ESLSKTWCAQ, from the internal region of the protein sequence according to NP_001010922.1.
Peptide Sequence	ESLSKTWCAQ
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	440603
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Preliminary experiments in Human Ileum and Testis lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

