

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-BCL2L15 ANTIBODY**

**SKU:** EB10383

## **SPECIFICATIONS**

**Unit Size** 100 μg

Synonym /

Alias BCL2L15|pro-apoptotic Bcl-2 protein|FL|22588|C1orf178|Bfk|BCL2-like 15

**Names** 

**Accession** NP\_001010922.1 ID

**Blocking** EBP10383

**Peptide** 

Peptide with sequence ESLSKTWCAQ, from the internal region of the protein sequence according to 

**Peptide** 

**ESLSKTWCAQ** Sequence

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Human

440603 **Gene ID** 

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring\_medium.png

Grade **ELISA** 

Detection Antibody detection limit dilution 1:1000.

Limit

**Blot** 

Western Preliminary experiments in Human Ileum and Testis lysates gave no specific signal but low background (at

antibody concentration up to 1µg/ml).

**Application** 

Pep-ELISA Type