


GOAT ANTI-BCL2L15 ANTIBODY

SKU: EB10384



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names BCL2L15|pro-apoptotic Bcl-2 protein|FLJ22588|C1orf178|Bfk|BCL2-like 15

Accession ID NP_001010922.1

Blocking Peptide EBP10384

Immunogen Peptide with sequence C-SPTLQVASRN, from the internal region of the protein sequence according to NP_001010922.1.

Peptide Sequence C-SPTLQVASRN

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human
Reactive Species	Human
Human Gene ID	440603
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Brain (Cerebellum). Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

