

# GOAT ANTI-BDH2 / DHRS6 (AA 92 TO 103) ANTIBODY

SKU: EB08340



---

## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias Names</b>	oxidoreductase UCPA dehydrogenase/reductase (SDR family) member 6 UNQ6308 UCPA-OR PRO20933 FLJ13261 EFA6R DHRS6 3-hydroxybutyrate dehydrogenase, type 2 BDH2
<b>Accession ID</b>	NP_064524.3
<b>Blocking Peptide</b>	EBP08340
<b>Immunogen</b>	Peptide with sequence DCEEKDWD FSMN, from the internal region of the protein sequence according to NP_064524.3.
<b>Peptide Sequence</b>	DCEEKDWD FSMN
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human, Mouse, Rat, Cow
<b>Human Gene ID</b>	56898
<b>Mouse Gene ID</b>	69772
<b>Rat Gene ID</b>	295458
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:128000.
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS

- [Data Sheet](#)