

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-BDH2 / DHRS6 (AA 92 TO 103) ANTIBODY

SKU: EB08340

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

oxidoreductase UCPA|dehydrogenase/reductase (SDR family) member 6|UNQ6308|UCPA-**Alias** OR|PRO20933|FLJ13261|EFA6R|DHRS6|3-hydroxybutyrate dehydrogenase, type 2|BDH2 **Names**

Accession

NP 064524.3 ID

Blocking

EBP08340 **Peptide**

Peptide with sequence DCEEKDWDFSMN, from the internal region of the protein sequence according to Immunogen

NP 064524.3.

Peptide

DCEEKDWDFSMN Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Cow **Species**

Human

56898

Gene ID Mouse

69772

Gene ID Rat Gene ID 295458

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:128000.

Limit

Application Pep-ELISA

Type