

GOAT ANTI-BHLHE22 ANTIBODY

SKU: EB12109



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	BHLHE22 basic helix-loop-helix family, member e22 BHLHB5 Beta3 CAGL85 TNRC20 basic helix-loop-helix domain containing, class B, 5 class B basic helix-loop-helix protein 5 class E basic helix-loop-helix protein 22 trinucleotide repeat containing 20 trinucleotide repeat-containing gene 20 protein
Accession ID	NP_689627.1
Blocking Peptide	EBP12109
Immunogen	Peptide with sequence C-HKSLSASTSKRLE, from the internal region (near N Terminus) of the protein sequence according to NP_689627.1.
Peptide Sequence	C-HKSLSASTSKRLE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	27319
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 37kDa band observed in Human (Cerebellum) Brain lysates (calculated MW of 37kDa according to NP_689627.1). Recommended concentration:0.1-0.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

