

GOAT ANTI-BICC1 ANTIBODY

SKU: EB09817



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	BICC1 DKFZp686I21119 bicaudal C homolog 1 (Drosophila) bicaudal C homolog 1
Accession ID	NP_001073981.1
Blocking Peptide	EBP09817
Immunogen	Peptide with sequence C-QHISQTYNISVS, from the internal region of the protein sequence according to NP_001073981.1.
Peptide Sequence	C-QHISQTYNISVS
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Human Gene ID	80114
Mouse Gene ID	83675
Rat Gene ID	361832
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)

