

GOAT ANTI-BIN3 ANTIBODY

SKU: EB08236



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	MEPCE[methylphosphate capping enzyme]BCDIN3[FLJ20257]bin3, bicoid-interacting 3 bin3, bicoid-interacting 3, homolog
Accession ID	NP_062552.2
Blocking Peptide	EBP08236
Immunogen	Peptide with sequence C-TPHNTSKGFQRP, from the internal region of the protein sequence according to NP_062552.2.
Peptide Sequence	C-TPHNTSKGFQRP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	56257
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx. 75kDa band observed in lysates of cell line K562 (calculated MW of 74.4kDa according to NP_062552.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

