



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-BIN3 ANTIBODY

SKU: EB08236



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa



Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

MEPCE|methylphosphate capping enzyme|BCDIN3|FLJ20257|bin3, bicoid-interacting 3|bin3, bicoid-interacting

Alias 3, homolog Names

Accession

NP_062552.2

Blocking

EBP08236 **Peptide**

Peptide with sequence C-TPHNTSKGFQRP, from the internal region of the protein sequence according to **Immunogen**

NP 062552.2.

Peptide

C-TPHNTSKGFQRP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

56257 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Western Approx. 75kDa band observed in lysates of cell line K562 (calculated MW of 74.4kDa according to

Blot NP_062552.2). Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB Type

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa